



Reference : 06-607CASE

Scharlau Microbiology - Technical data sheet

Product :  
BASIC FUCHSIN (250) SELECTIVE SUPPLEMENT



F

### Specificacion

Selective and differential supplement for the coliform micro organisms detection.

### Packaging content

The box contains 10 vials. Each vial is sufficient to supplement 500 mL of Endo Agar Base Ref. 01-589.

### Vial content

Necessary amount for 500 mL of solid medium.

Basic fuchsine.....250 mg

Ethanol.....5 mL

### Directions

Mix the liquid with the powder by pressing down on the cap. Shake to dissolve and aseptically add the solution to the boiled medium. Homogenize and use as per each monograph medium.

**Note:** Do not re-melt the solid media.

### Precautions

-This product should be for laboratory use only.

-Does not use beyond stated expiry date.

### Applicable media

Ref. 01-589 Endo Agar Base

### Sterility control

No growth to 7 days at 37°C into Thioglichollate USP and TSB.

### Storage

For laboratory use only. Keep tightly closed, away from bright light, in a cool dry place (+2°C to 8°C).