



Specification

Selective solid medium used for the detection of *Pseudomonas aeruginosa* according to the EN 12780-2002 and ISO 16266 Standard.

Presentation

	Packaging Details	Shelf Life	Storage
30 Prepared Plates 55 mm Plates for filtration purposes with: 9 ± 1 ml	1 box containing 3 plastic bags with 3 RD-PACK. Each package contains 10 plates of 55 mm.	6 months	2-25°C

Composition

Composition (g/l):	
Gelatin Peptone.....	16.000
Casein Peptone.....	10.000
Magnesium chloride.....	1.400
Potassium sulfate.....	10.000
Cetrimide.....	0.200
Sodium Nalidixate.....	0.015
Glycerol.....	10.000 ml
Agar.....	13.600

Description /Technique

Description:

The CN Selective Medium for *Pseudomonas* was progressively developed from the basic media of King, Ward and Raney for the production of pigments. Browne and Lowbury add the cetrimide as selective agent and Goto and Enomoto improves efficiency by adding nalidixic acid. The presence of both inhibitors eliminates the contaminant microbiota from heavily polluted specimens and was adopted by the CEN (Centre Europe ende Normalisation) in its EN Standard 12780 for the detection of *Ps. aeruginosa* by filtering membrane in water.

Technique:

Collect, dilute and prepare samples and volumes to be filtered as required according to specifications, directives, official standard regulations and/or expected results.

Filter the sample through a 0.45 mm Ø pore membrane and apply it onto the surface of the agar.

Incubate the plates air room atmosphere at $35 \pm 2^\circ\text{C}$ for 24-48h.

(Incubation times longer than those mentioned above or different incubation temperatures may be required depending on the sample, on the specifications).

After incubation, count the colonies with a blue-greenish appearance due to pigment production by *Pseudomonas* sp.

Calculate total microbial count per ml of sample by multiplying the average number of colonies per plate by the inverse dilution factor.

Report results as Colony Forming Unit (CFU) per ml along with incubation time and temperature.

Presumptive isolation of *Pseudomonas* sp must be confirmed by further tests.

Quality control

Physical/Chemical control

Color : Off-white / opalescent pH: 7.2 ± 0.2 at 25°C

Microbiological control

Membrane Filtration /Practical range 100 ± 20 CFU; Min. 50 CFU (Productivity)./ 10^4 - 10^6 CFU for Selectivity.

Microbiological control according to ISO 11133:2014

Aerobiosis. Incubation at $36 \pm 2^\circ\text{C}$, reading at 44 ± 4 h

Microorganism

Ps. aeruginosa ATCC® 9027, WDCM 00026

Ps. aeruginosa ATCC® 27853, WDCM 00025

Ps. aeruginosa ATCC® 10145, WDCM 00024

Escherichia coli ATCC® 8739, WDCM 00012

Enterococcus faecalis ATCC® 29212, WDCM 00087

Growth

Good (≥ 50 %)

Good (≥ 50 %)

Good (≥ 50 %)

Inhibited

Inhibited

Sterility Control

Incubation 48 hours at 30 - 35°C and 48 hours at 20 - 25°C : NO GROWTH

Check at 7 days after incubation in same conditions

Bibliography

- BROWN, V.L. & E.J.L. LOWBURY (1965) Use of an improved Cetrimide Agar Medium and of culture methods for *P. aeruginosa*. J., Clin. Pathol. 18:752.
- GOTO S. & S. ENOMOTO (1970) Nalidixic acid cetrimide agar. A new selective plating medium for the selective isolation of *P. aeruginosa*. Jpn. J. Microbiol. 14:65.
- ISO 11133:2014. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- ISO 12780:2002 Water Quality. Detection and enumeration of *Pseudomonas aeruginosa* by membrane filtration.
- ISO 16266 Standard (2006) Water Quality. - Detection and enumeration of *Pseudomonas aeruginosa*. - Method by membrane filtration.
- KING, E.O., M.K. WARD & E.E. RANEY (1954) Two simple media for the demonstration of pyocyanin and fluorescein. J. Lab. Clin. Med. 44:301.
- ROBIN, T. & J.M. JANDA (1984) Enhanced recovery of *P. aeruginosa* from diverse clinical specimens on a new selective agar. Diag. Microbiol. Infect Dis. 2:207.
- SCHWEIZERISCHE LEBENMITTELSBUCH (2005) Kap. 56 Mikrobiologie. Bundesamt für Gesundheit. Direktionsbereich Verbraucherschutz. Bern.