

**Iron(II) ethylenediammonium sulfate tetrahydrate, secondary standard for titrations, Titrasure®****Identification**

$\text{FeC}_2\text{H}_4(\text{NH}_3)_2(\text{SO}_4)_2 \cdot 4\text{H}_2\text{O}$

M = 382,14 g/mol

CAS [113193-60-5]

Taric code: 2921 21 00

Specifications

assay (on dried sample).....min. 99,5 %

This volumetric standard is directly traceable to NIST Sodium Oxalate (Reductometric Standard).

Physical data

• Bulk density: ~ 600 kg/m³

Transport/storage

• 10°C - 30°C